

Command	ADESTROY
PURPOSE	Destroy analyzer and free the occupied space
PARAMETERS	
ANLID	Analyzer identifier
/NOCONFIRM	No confirmation of analyzer destruction
/NOLIST	Suppress output of command
FUNCTION	The specified analyzer is destroyed and the space it occupies is freed. To reset spectra contents to zero use ACLEAR!
REMARKS	A specific analyzer or a group of analyzers can be selected by use of the wildcard “*”. Resident analyzers (analyzers defined in the analysis program by the macro AGEN) cannot be destroyed. If not suppressed by /NOCONFIRM, prompting for confirmation precedes any analyzer destruction in the case of wildcards used.
EXAMPLE	ADES * Destroy all dynamically created analyzers.